



Schema Changes in ICAT 6.2 and 7.0

Hints for a seamless upgrade

Rolf Krahl 

ICAT F2F Meeting, 10 February 2026, Abingdon, UK

Schema Changes in ICAT 6.2

- Schema Changes in ICAT 6.2 (released 2025-11-28):
 - Add a new entity type Subject to add keywords to data publications #327
 - Add startDate and endDate fields to Instrument #349
 - Add internalId attribute to DataPublication #355
 - Add acknowledgement attribute to FundingReference #369
 - Increase size of location column for Datafile and Dataset #364
- Only incremental changes, no compatibility issues.
- Upgrade: recommend to run the SQL upgrade script to keep the size of location columns in the database consistent with the schema.
- Besides that: just start to use the new entity type and new attributes if you need them.

Migrating from DataPublication.subject to Subject

- Assuming you store a semicolon separated list in DataPublication.subject, you may migrate to the new Subject entity:

```
for pub in client.search("DataPublication"):  
    for subj in pub.subject.split(';'):  
        subject = client.new("Subject", dataPublication=pub)  
        subject.name = subj.strip()  
        subject.create()
```

- You still need to update your scripts to use DataPublication.subjects in place of DataPublication.subject
- You will need appropriate access permissions for class Subject.
- You may now use the new attributes in class Subject.

Schema Changes in ICAT 7.0

- Schema Changes in ICAT 7.0 (planned):
 - Make the relationship between `Sample` and `Investigation` many-to-many #231
 - Add a `pid` property to `SampleType` and change the uniqueness constraint #326
- Incompatible changes. You must run the SQL upgrade script!
- It is recommended to prepare the upgrade with respect to the `Sample` table while the old ICAT server is still running, see subsequent slides.

Changes in SampleType

- Changes to SampleType are:
 - drop the relation to Facility,
 - add new pid and description attributes,
 - allow NULL in molecularFormula,
 - change the uniqueness constraint to be just pid.
- Caveat for the upgrade: the new pid attribute must be populated with unique values. The upgrade script cares about that, setting dummy values.
- If you don't care about PIDs for sample types, you can just keep using dummy values and disregard them.

Changes in SampleType (cont.)

- You may need to update your scripts, mostly for creating new SampleType:

```
sampletype = client.new("SampleType")
if 'pid' in sampletype.InstAttr:
    # ICAT schema >= 7.0
    sampletype.pid = str(uuid.uuid4())
sampletype.name = "Some_SampleType_name"
sampletype.molecularFormula = "Some_formula"
if 'facility' in sampletype.InstRel:
    # ICAT schema < 7.0
    sampletype.facility = client.assertedSearch("Facility")[0]
sampletype.create()
```

Changes in Sample

- Changes to Sample are:
 - drop the many-to-one relation from Sample to Investigation,
 - add class InvestigationSample to establish a many-to-many relation from Sample to Investigation,
 - change the uniqueness constraint in Sample to be just pid.
- Caveat for the upgrade: the pid attribute must be populated with unique values.
- Additional challenge here: the pid attribute did exist in the old schema. Depending on the existing values, there may be obstacles for the upgrade: non-unique values, values that may conflict with the dummy values that the upgrade script tries to set.
- Recommended strategy: prepare the upgrade in the running ICAT, populate the pid attribute before trying to upgrade.
- We provide a script helping with that, see next slide.

Preparing the Upgrade

- Helper script `check-samples`, see `icatproject-contrib/scripts#1`, implementing four subcommands:
 - `stats` Display some statistics about the values in `Sample.pid` and raises warnings if there are any obstacles for the upgrade,
 - `lsdup` List all duplicate `Sample.pid` values,
 - `setpids` populate `Sample.pid`, e.g. set unique dummy values for all `Sample` where `pid` is `NULL`,
 - `dedup` deduplicate `Sample.pid` values, e.g. change duplicate values so that they are unique by appending a suffix.
- If `check-samples stats` does not issue a warning, the upgrade via the SQL script should work.
- Still, it is recommended to run `check-samples setpids` if there are any `NULL` values, just to be safe.

Changes in Sample (cont.)

- You may need to update your scripts, for instance creating a new Sample related to a given Investigation:

```
sample = client.new("Sample")
sample.name = "Some_name"
sample.pid = str(uuid.uuid4())
if 'investigationSamples' in sample.InstMRel:
    # ICAT schema >= 7.0
    inv_samp = client.new("InvestigationSample",
                          investigation=investigation)
    sample.investigationSamples.append(inv_samp)
else:
    # ICAT schema < 7.0
    sample.investigation = investigation
sample.type = samplotype
sample.create()
```

- You will need appropriate access permissions for class InvestigationSample.