



Getting started – how to set up the common data tables

NOBUGS 2012 ICAT Workshop
27th September 2012

Kevin Phipps
Scientific Computing Department, RAL

What needs to be in place?

- ICAT (4.2) installed and running on an application server such as Glassfish
- An authentication plugin
 - 2 standard versions available: DB and LDAP
 - Download from <http://www.icatproject.org/mvn/site/>



Using the database authentication plugin

- When the authentication ear file is deployed a database table `PASSWD` is created
- 2 columns `USERNAME` and `ENCODEDPASSWORD`
- Add a username and a password – the password must be in plain text and not encoded as the column name suggests. This will change in a later version.



Connecting to ICAT

- ICAT runs as a web service and is described by its WSDL document such as the one at:
<http://www.icatproject.org:8080/ICATService/ICAT?wsdl>
- To connect using Java use the client library at:
<http://www.icatproject.org/mvn/site/icat/4.2.0/icat.client/installation.html>
- Alternatively generate your own client using a tool such as NetBeans or connect in Python via the suds library by pointing it at the WSDL URL



Making the initial connection

```
URL hostUrl = new URL("http://localhost:8080");
URL icatUrl = new URL(hostUrl, "/ICATService/ICAT?wsdl");
QName qName = new QName("http://icatproject.org", "ICATService");
ICATService service = new ICATService(icatUrl, qName);
ICAT icat = service.getICATPort();
```

- Keep hold of the reference to “icat”. You will need it for every call you want to make to the server.



Setting up login credentials

```
Credentials credentials = new Credentials();  
List<Entry> entries = credentials.getEntry();  
Entry e;  
e = new Entry();  
e.setKey("username");  
e.setValue("admin");  
entries.add(e);  
e = new Entry();  
e.setKey("password");  
e.setValue("secret");  
entries.add(e);
```



Login

- Pass the credentials to the login method specifying the authentication mechanism to use

```
String sessionId = icat.login("db", credentials);
```

- Hold on to the sessionId. It is required when making subsequent calls to ICAT.



Setting up basic authorisation

- Add the admin user you have set up to ICAT

```
User adminUser = new User();  
adminUser.setName("admin");  
adminUser.setId(icat.create(sessionId, adminUser));
```

- Set up an ICAT Group

```
Group rootGroup = new Group();  
rootGroup.setName("root");  
rootGroup.setId(icat.create(sessionId, rootGroup));
```

- Add the admin user to that group

```
UserGroup userGroup = new UserGroup();  
userGroup.setUser(adminUser);  
userGroup.setGroup(rootGroup);  
icat.create(sessionId, userGroup);
```



Assigning privileges


- Create a rule – the example below gives our user “admin” Create, Read, Update and Delete privileges on Facility objects

```
Rule rule = new Rule();  
rule.setGroup(rootGroup);  
rule.setWhat("Facility");  
rule.setCrudFlags("CRUD");  
icat.create(sessionId, rule);
```

- Add further rules for each type of object to be created




Schema - The main components



- **Facility – an experimental facility**
ICAT must contain one before you can add any data.
- **Investigation – an investigation or experiment**
This is typically defined by a group of users visiting an experimental facility for a defined period of time during which they create data they want to save.



Main components (cont.)



- **Dataset** – a collection of data files and parameters
Typically contain a number of DatasetParameters (to record settings, times, temperatures etc) and a list of associated Datafiles.
- **Datafile** – a file and associated parameters
Examples of these are configuration files, image files, binary data files. Each Datafile can have a list of DatafileParameters associated with it.



Creating a Facility in ICAT

```
Facility facility = new Facility();  
facility.setName("ISIS"); // required  
facility.setFullName("ISIS Pulsed Neutron and Muon Source"); // optional  
facility.setId(icat.create(sessionId, facility));
```



Creating an Investigation

- An InvestigationType must be created first

```
InvestigationType invType = new InvestigationType();  
invType.setFacility(facility);  
invType.setName("calibration");  
invType.setDescription("A calibration experiment");  
invType.setId(icat.create(sessionId, invType));
```

- Then use it to create the Investigation

```
Investigation inv = new Investigation();  
inv.setFacility(facility);  
inv.setName("Investigation 1");  
inv.setTitle("A more grand title for Investigation 1");  
inv.setType(invType);  
inv.setId(icat.create(sessionId, inv));
```



Creating a Dataset

- A DatasetType must be created first

```
DatasetType rawDsType = new DatasetType();  
rawDsType.setFacility(facility);  
rawDsType.setName("experiment_raw");  
rawDsType.setDescription("RAW data collected during an experiment");  
rawDsType.setId(ocat.create(sessionId, rawDsType));
```

- Then use it to create Datasets of that type
(Covered in more detail later in the Workshop)



Creating a Datafile

- A DatafileFormat must be created first

```
DatafileFormat nexusFormat = new DatafileFormat();
nexusFormat.setFacility(facility);           // required
nexusFormat.setName("nexus");               // required
nexusFormat.setVersion("4.3.0");           // required (I think!)
nexusFormat.setType("HDF5");
nexusFormat.setDescription("ISIS Muon format");
nexusFormat.setId(this.icat.create(this.sessionId, nexusFormat));
```

- Then use it to create Datafiles of that type
(Covered in more detail later in the Workshop)



Adding parameters to objects

- A `ParameterType` must be created first

```
ParameterType protonChargeType = new ParameterType();  
protonChargeType.setFacility(facility); // required  
protonChargeType.setName("total_proton_charge"); // required  
protonChargeType.setUnits("uAh"); // required  
protonChargeType.setValueType(ParameterValueType.NUMERIC); // required  
protonChargeType.setApplicableToDatafile(true);  
protonChargeType.setApplicableToDataset(true);  
protonChargeType.setApplicableToInvestigation(true);  
protonChargeType.setApplicableToSample(true);  
protonChargeType.setDescription("The total proton charge on target during the  
    collection period");  
protonChargeType.setUnitsFullName("micro amp hours");  
protonChargeType.setEnforced(true);  
protonChargeType.setMinimumNumericValue(new Double(0));  
protonChargeType.setVerified(true);  
protonChargeType.setId(ocat.create(sessionId, protonChargeType));
```



Basic setup complete

- ICAT contains a Facility, an Investigation, some DatasetTypes, ParameterTypes and DatafileFormats
- ICAT is now ready to ingest Datasets and Datafiles
- Further configuration required for a more extensive setup using FacilityCycles, Instruments, Applications etc but they are all done in the same way



And finally...

- Any questions?
- Discussion about possibility of agreeing a basic set of Types and Formats?

- And thanks - to Tom Griffin for providing the ISIS related examples used in this presentation



Dataset_Type

NAME	DESCRIPTION
analyzed	Analyzed data
experiment_cal	RAW data collected at the facility during a calibration run.
experiment_eng	RAW data collected at the facility during an engineering test.
experiment_raw	RAW data collected at the facility during an experiment.
reduced	Reduced Data
simulation	Simulation data
special_cal	"Calibration data not acquired through the normal Data acquisition system.

Investigation_Type

NAME	DESCRIPTION
calibration	A set of Calibration
commercial_experiment	A scientific experiment performed by a commercial company
engineering	"Calibration, first light data, alignment, ?"
experiment	A scientific experiment.
simulation	

Datafile_Format

NAME	VERSION	FORMAT_TYPE	DESCRIPTION
nexus	3.0.0	HDF5	Neutron and X-Ray data format.
isis neutron raw v8	8	binary	
isis muon raw v4	4	binary	
nexus	1	binary	
nexus	2.1.0	HDF4	ISIS Muon format
isis neutron raw	2	binary	
nexus	4.3.0	HDF4	ISIS Muon format
nexus	4.1.0	HDF4	ISIS Muon format

